

Issue 4 – 02.03.2020 B-Lab-For-009a

Page 1 of 2

BioLabTests

Unit 3 Parade Court Central Boulevard, Prologis Park Coventry, CV6 4QL United Kingdom

Tel: +44 (0) 333 240 8308 Fax: +44 (0) 2476 338081 info@biolabtests.com www.biolabtests.com

Certificate of Analysis

Customer Name: PG BISON Date Received: 20.07.2020
Customer Contact: Wouter Britz Date Analysed: 28.07.2020
Customer Email/Phone: wbritz@pgbison.co.za Date Reported: 24.08.2020

Certificate Number: BL050/2020

Tick box if customer has provided untreated control samples for inclusion in testing

Test Method: Determination of Antibacterial Activity using Test Based on MOD ISO 22196: 2011

BioLabTests Sample		Contact Time*		Reduction against Initial		
Reference Number	Sample	Test Bacteria	0 Hrs	24 Hrs	Log ₁₀	%
BL/039	MELAMINE-FACED BOARD WITH ANTIMICROBIAL ADDITIVE	MRSA	8.13 x10 ⁴	≤100	≥2.91	≥99.88%
		E. coli	7.81 x10 ⁴	≤100	≥2.89	≥99.87%

^{*}Numbers represent Colony Forming Units at representative contact times.

Data show the changes in size of the initial MRSA (ATCC® 43300™) and *E. coli* (ATCC®8739™) populations due to their contact with the treated surface of the samples. The samples detailed in this CoA were tested for 24 hours (±1) at 35°C (±2) under a relative humidity of >90% as specified in the ISO 22196 method. The results detailed in this CoA relate only to the items tested. This certificate shall not be reproduced except in full, without written approval of BioLabTests.

Please note; sample was received in a ready to test state from the customer apart from sterilisation which is noted below. All testing is performed on site at the BioLabTests address above unless otherwise disclosed.





Page 2 of 2

BioLabTests

Unit 3 Parade Court Central Boulevard, Prologis Park Coventry, CV6 4QL United Kingdom

Tel: +44 (0) 333 240 8308 Fax: +44 (0) 2476 338081 info@biolabtests.com www.biolabtests.com

Noted during testing				
	Note			
No. of samples supplied per test:	6 of 6 supplied			
Sample size:	Average size: 50 x 50mm			
Sample sterilisation method:	DRY HEAT AT 65°C FOR ONE HOUR			
Film size:	Average size: 40 x 40mm			
Inoculation volume:	Volume: 0.1ml			
	MRSA: 8.13 x10 ⁵ CFU/ml			
Initial bacterial concentration:	<i>E. coli</i> : 7.81 x10 ⁵ CFU/ml			

The material of the cover slips is polyethylene with a thickness of 0.045mm.

Please note – For all information regarding deviations from ISO 22196:2011, please request statement B-Lab-State-007

Ria Warren BSc (Hons), MRes Quality Manager, BioLabTests Ltd

